



Marine Mammal Monitoring System



- This system has been using for Monitoring and record the sound of marine mammals. Sound data can be saved to SD card, and simultaneous transmitting to the base station using wireless LAN. It can also get more accurate positioning data by using two or more tracking units.

SPECIFICATION

Hydrophone Array	4 Hydrophones (2 sets of horizontal arrays)
Span of Array	200mm(horizontal), 100mm(vertical)
Receiving Frequency Range	10kHz~200kHz
Battery Type	Lithium Battery
Maximum Detectable Range	100m
Outline Dimensions	
Hydrophone Array	φ 224mm×L700mm
Electronic Device	φ 132mm×L460mm
Weight in Air	6kg